Engineered for you

Orbital Shaker MO 10

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info@hettichbenelux.com

Content

1. Basic Data.....	4
Safety.....	4
Warranty	5
Disclaimer of Warranty	5
Important Information	5
Environmental conditions	5
2. First Steps.....	6
Setup.....	6
3. Operation.....	6
Switching on	6
Parameter setting.....	7
Operation with timer.....	7
Automatic re-start after electrical failure.....	8
4. Assembly of Accessories	8
Platforms	8
5. Troubleshooting	8
6. Maintenance	9
Cleaning	9
Service	10
7. Technical Data	11
8. Article description accessories	12
9. Explanations	13

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1. Basic Data

The HLC by Hettich Benelux Orbital Shaker MO 10 is used for the mixing of solutions in closed reaction tubes, plates and flasks up to 10 kg. If there are any uncertainties regarding the liquids, which should be shaken, please contact Hettich Benelux. The HLC by Hettich Benelux Orbital Shaker can be integrated in an incubator and therefore is perfectly suited for applications, where test samples must be tempered. The HLC by Hettich Benelux instrument is intended for indoor usage only. Please use only HLC by Hettich Benelux accessories or accessories recommended by HLC by Hettich Benelux.

Safety

The instrument complies with the standards and directives mentioned in the applicable CE declaration.

Please take note of the following safety measures:

- » Do only connect the instrument to an earthed mains power socket of 115 V - 230 V, 50 – 60 Hz.
- » The objects to shake must be put on the tray in the most centred position possible
- » The installation surface of the instrument must be as even as possible and may not be slippery so that the instrument can fulfil the standard IN 61010-2-051
- » If liquid gets into the instrument, immediately pull out the plug and contact our service department to ensure complete safety.
- » Contact with highly flammable fluids must be avoided.
- » The ventilation slots of the device may not be blocked to ensure proper ventilation at any time.
- » If the mains cable is damaged, it should only be replaced by an equal one.

Warranty

All functions of the instrument have been tested thoroughly. After that, the instrument and the accessories were packed carefully in perfect condition.

If, however, any damages or defects are detected during installation or setup, please inform Hettich Benelux immediately (see our Terms and Conditions at www.hettichbenelux.com). If used under normal laboratory conditions according to the Operating Manual, Hettich Benelux grants a warranty of two years for material and workmanship, starting with the day of shipment.

Disclaimer of warranty

In case of unintended use of the device, the manufacturer assumes no responsibility.

Important Information

This manual complies with the standards and directives mentioned in the EN DIN EN 61010-1 standard. IQ, OQ and PQ documents are available on request.

If the Orbital Shaker is operated with a high shaking frequency, the vibrations of the instrument might transmit to the under-ground. Please do consider this when placing the instrument.

Environmental conditions

The ideal operation environment of the is 80% relative air humidity or below and between +0°C - +65°C ambient air temperature. These values should neither be underrun nor exceeded.

2. First Steps

Setup

Before inserting the plug, take care that the voltage of the plug corresponds to the voltage of the instrument. The power socket used must be earthed (protective earth conductor contact). If these procedures are not followed, safe operation cannot be guaranteed and/or the instrument may be damaged. Observe the ambient conditions (temperature, humidity, etc.) listed in Technical Data.

Plug the mains cable into the IEC socket on the back of the instrument and connect it with the 115V-230 V, 50-60 Hz via an earthed mains power socket.

The Orbital Shaker has 4 feet for regulating the height in order to even out the surface as much as possible.

3. Operation

Switching on

For shaking, the mains switch on the right side on the front of the instrument must be switched on.

When activating the mains switch, the LCD display will light up, the instrument will perform an auto check and the display will show the version of installed software.

Once the instrument is initialized, the display will show the following:

First row: shows the actual speed and the elapsed time (initially zero)

Second row: shows the programmed speed and time of the last session. In the middle of the row, it is indicated if the instrument is activated (ON) or not (##).

Parameter setting

The cursor on the display initially is located in selected RPM. Press **MIN/RPM** to change the cursor from RPM to MIN (time).

By means of the selection switches (**▲** and **▼** buttons), adjust the desired values. Note: if you push the **▲** and **▼** buttons simultaneously when the cursor is in RPM position, the instrument will be adjusted to its maximum available speed.

Touch the **START/STOP** key and the display will show (ON), the instrument will initiate shaking in a soft and gradual way.

To stop shaking, touch **START/STOP** again. The display will now show (##).

Note: you can adjust the speed of the instrument while it is working. For doing so, use the **▲** and **▼** buttons.

***Safety system** for motor protection: if the instrument surpasses maximum motor power, it will be stopped and the message **STOP** is blinking on the display. Press the **START/STOP** button in order to stop the instrument correctly, and once you have solved the problem that causes these adverse conditions, you can start the instrument again.

Operation with timer

First of all, adjust the desired value of RPM with the selection switches (**▲** and **▼** buttons).

Press the **MIN/RPM** button, and the cursor will be located in the selected time. By means of the selection switches (^ and v buttons) adjust the desired minutes. Note: if you select "0" minutes, the instrument will operate continuously. Press the **START/STOP** button, and the display will show (ON), and it will start shaking softly and gradually, increasing the speed until it reaches the selected one.

If the programmed minutes are reached, the instrument will stop automatically.

Automatic re-start after electrical failure

In case of accidental cut off of the electrical supply, the instrument returns to "ON" with the same adjusted parameters when the power returns. In the display it will show the symbol "!" to indicate this situation.

Note: disconnect the instrument correctly; push **START/STOP** button and then the mains switch. If it is not done this way, the instrument will understand that an electrical failure has happened.

4. Assembly of Accessories

Platforms

The platforms can be placed on the top of the instrument. Use the black safety screws for fixing the platform. Check that the screws are tight. For removing the platform, loosen the screws again and pull up the platform. Never use the instrument without checking if the platform is fixed to the shaker.

5. Troubleshooting

The display remains blank

Please check if the mains switch on the front side is switched on. If it is switched on but no voltage is present at the socket, check the micro-fuse and

replace if required. (IEC 127-2/III, 250 V, 2 A time-lag). This fuse - and a spare fuse - are located in the IEC-bushing (in which the connection cable is plugged). The fuse box can be pulled out with a screwdriver.

6. Maintenance

Cleaning

Regularly clean the housing and the platforms of the Orbital Shaker.

Precautions for avoiding electric shock

Electronic devices can cause electric shocks in case of an operating error. Never try to repair electric parts. Never open the housing.

- Switch off the instrument and disconnect it from the power supply before starting with cleaning or disinfection works.
- Never let get liquids inside the housing (ventilation slit).
- Do not perform spray disinfection.
- Do only connect the instrument with the power supply if it is completely dry.

The repair service may only be performed by authorized staff trained by the manufacturer. A modification of the instrument is not permitted.

Caution when handling aggressive chemicals

Do not use aggressive chemicals like e.g. strong and weak bases, strong acids, formaldehyde, acetone, halogenated hydrocarbons or phenol for cleaning the instrument and its accessories.

- In case of contamination with aggressive chemicals, clean the instrument with a neutral detergent immediately.
- Use neither corrosive detergents nor aggressive solvents or abrasive polishing agents.

Cleaning

1. Please disconnect the instrument from the power supply before you start cleaning.
2. Please clean all outer parts of the instrument with a mild detergent and a lint-free cloth.
3. Please wipe off the detergent with distilled water.
4. Please dry all cleaned parts.

Disinfection

1. Please disconnect the instrument from the power supply before you start disinfecting.
2. Please let the instrument cool down.
3. Please clean the instrument as described above.
4. Select a disinfection method compliant to the applicable local legal regulations and directives.
5. Please wipe off all outer parts of the instrument with the disinfectant and a lint-free cloth.

Service

If a technical problem arises, please contact your local dealer or the Hettich Benelux service department. The Hettich Benelux contact information can be found at www.hettichbenelux.com. If required, the instrument will be returned for repair – please take note of our service guidelines, which can be found at www.ditabis.com as well.

Decontamination before shipment

If you send the instrument to the authorised technical service for repair or to your distributor for disposal, decontaminate all parts you want to send. Document the decontamination in a Decontamination Certificate (incl. serial number) and include it with the shipment.

7. Technical Data


Technical Data	Orbital Shaker MO 10
Shaking frequency	50 - 400 rpm
Orbit	10 mm round
Max. Load	10 kg
Timer	0 – 120 min
Dimensions W x D x H	410 x 570 x 140 mm
Platform size	350 x 400 mm
Weight	13,5 kg
Power input	50 W
Electr. supply	110...230 V, 50...60 Hz
Fuse	F1A
Protection class	IP21
Environmental conditions	
Ambient air temperature in service (non-condensing)	0° – 65°C
Relative humidity	max. 80%

8. Article description accessories

Platforms and Accessories

Art. No.	Article
300200025	Non-slipping Platform, 350 x 400 mm
300200026	Double Non-slipping Platform, 350 x 400 mm
300200027	Universal Platform with 5 rollers, 350 x 400 mm
300200028	Double Universal Platform with 10 rollers , 350 x 400 mm
300200029	Roller for Universal Platform
300200030	Platform fastening for Erlenmeyer clips
300200031	Erlenmeyer Clip 50 ml, number per platform 40 x 50 ml Erlenmeyer Clips
300200032	Erlenmeyer Clip 100 ml, number per platform 20 x 100 ml Erlenmeyer Clips
300200033	Erlenmeyer Clip 250 ml, number per platform 12 x 250 ml Erlenmeyer Clips
300200034	Erlenmeyer Clip 500 ml, number per platform 9 x 500 ml Erlenmeyer Clips
300200035	Erlenmeyer Clip 1000 ml, number per platform 9 x 1000 ml Erlenmeyer Clips
300200036	Erlenmeyer Clip 2000 ml, number per platform 5x 2000 ml Erlenmeyer Clips
300200037	Platform for 5 separating funnels 250 ml, 500 ml; 3 x 1000 ml
300200038	Platform for 6 separating funnels 250 ml, 500 ml

9. Explanations

	Caution! Risk of electric shock!
	Caution!
	Important information
	Earth conductor
	Fuse

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