

ATLAS DUCT SPECIFICATIONS



For more information please visit:
www.labogene.com

Atlas Duct	Unit	Part no: 9.100.010.111
External dimensions (DxWxH)	mm	340 x 627 x 990
Working chamber, dimensions (DxWxH)	mm	285 x 573 x 535
Working opening (w x H)	mm	573 x 300
Exhaust airflow	m ³ /hr	300
Noise level, EN 11201	dB(A)	46
Light Intensity variable	Lux	700
HEPA filters, EN 1822 H-14 99.995% MPPS		Yes
Carbon filter		No
Voltage/frequency	V/Hz	220-240/50-60
Power consumption	amps	0,2
Fuses	A	10
Net weight	kg	25
Window material		Polycarbonate
Worktop		AISI 304 stainless steel
Cabinet material		Powder coated steel

LaboGene ApS
Industrivej 6-8
DK-3540 Lyngø

Tel (+45) 3940 2566
Fax (+45) 3995 2566
Mail info@labogene.com
Web www.labogene.com

