

# SIRIUS WITH 20 CM WORK OPENING

## SPECIFICATIONS



For more information please visit:  
[www.labogene.com](http://www.labogene.com)

Sirius with 20 cm work opening	Unit	Part no: 9.002.001.072
External dimensions (DxWxH)	mm	600 x 690 x 1365
Working chamber, dimensions (DxWxH)	mm	457 x 687 x 665
Working opening (WxH)	mm	670 x 665
Exhaust airflow	m <sup>3</sup> /hr	200
Downflow	m/s	0,4
Noise level, EN 11201	dB(A)	54
Light Intensity variable	Lux	700
HEPA filters, EN 1822 H-14 99.995% MPPS		Yes
Carbon filter		Option
Voltage/frequency	V/Hz	220-240/50-60
Power consumption	amps	0,2
Fuses	A	10
Net weight	kg	100
Window material		Polycarbonate
Worktop		Stainless steel 304
Cabinet material		Powder coated steel

LaboGene ApS  
 Industrivej 6-8  
 DK-3540 Lyngø

Tel (+45) 3940 2566  
 Fax (+45) 3995 2566  
 Mail [info@labogene.com](mailto:info@labogene.com)  
 Web [www.labogene.com](http://www.labogene.com)

