

SIRIUS WITH ARM-HOLES SPECIFICATIONS



For more information please visit:
www.labogene.com

Sirius with arm-holes	Unit	Part no: 9.002.000.072
External dimensions (DxWxH)	mm	600 x 690 x 1365
Working chamber, dimensions (DxWxH)	mm	457 x 687 x 665
Working opening (WxH)	mm	2 x 250mm armholes
Exhaust airflow	m ³ /hr	100
Downflow	m/s	0,4
Noise level, EN 11201	dB(A)	54
Light Intensity variable	Lux	700
HEPA filters, EN 1822 H-14 99.995% MPPS		Yes
Carbon filter		Option
Voltage/frequency	V/Hz	220-240/50-60
Power consumption	amps	0,2
Fuses	A	10
Net weight	kg	100
Window material		Polycarbonate
Worktop		Stainless steel 304
Cabinet material		Powder coated steel

LaboGene ApS
Industrivej 6-8
DK-3540 Lyngø

Tel (+45) 3940 2566
Fax (+45) 3995 2566
Mail info@labogene.com
Web www.labogene.com

