

TITAN PCR SPECIFICATIONS



For more information please visit:
www.labogene.com

Titan PCR	Unit	Part no: 9.100.010.011
External dimensions (DxWxH)	mm	340 x 627 x 990
Working chamber, dimensions (DxWxH)	mm	285 x 573 x 535
Working opening (w x H)	m	573 x 400
Air velocity, Downflow (Stepless adjustable)	m/s	0.28
Air velocity, deviation	+/- %	10
Noise level, EN 11201	dB(A)	49
Light Intensity variable	Lux	700
HEPA filters, EN 1822		99.999% efficient H14
UV		Yes
Voltage/frequency	V/Hz	220-240/50-60
Power consumption	amps	0,2
Fuses	A	10
Net weight	kg	25
Window material		Polycarbonate
Worktop		AISI 304 stainless steel
Cabinet material		Powder coated steel

LaboGene ApS
Industrivej 6-8
DK-3540 Lyngø

Tel (+45) 3940 2566
Fax (+45) 3995 2566
Mail info@labogene.com
Web www.labogene.com

